

The logo icon consists of a stylized 'R' shape formed by two concentric, rounded lines. The outer line is black and the inner line is blue, creating a blue-filled 'R' with a black outline.

**rehametrics**

quantified physical and cognitive rehabilitation

# total control during the session



pause / continue / stop the exercise

increase / decrease the difficulty

increase / decrease session time


enables control of each exercise without stopping the session

# keyboard commands



AVAILABLE COMMANDS		
Pause the session	ENTER	
At the "Pause" screen:		
Continue	RIGHT ARROW	
End session	LEFT ARROW	
Positioning screen	ENTER	
During the session:		
Increase exercise difficulty	UP ARROW	+1 level of difficulty
Decrease exercise difficulty	DOWN ARROW	-1 level of difficulty
Increase session time	RIGHT ARROW	+30 seconds
Decrease session time	LEFT ARROW	-30 seconds

# remote control commands



AVAILABLE COMMANDS		
Pause the session	OK	
At the "Pause" screen:		
Continue	RIGHT ARROW	
End session	LEFT ARROW	
Positioning screen	OK	
During the session:		
Increase exercise difficulty	UP ARROW	+1 level of difficulty
Decrease exercise difficulty	DOWN ARROW	-1 level of difficulty
Increase session time	RIGHT ARROW	+30 seconds
Decrease session time	LEFT ARROW	-30 seconds

Thank you



Microsoft Health  

---

Innovation Awards 2015  
WINNER

